PACKAGE INSERT

AlloWrap[®] DS AlloWrap[®] Dry

Double-Sided Human Amniotic Membrane

THIS ALLOGRAFT IS DERIVED FROM VOLUNTARILY DONATED HUMAN TISSUES. IT IS INTENDED FOR SINGLE PATIENT, SINGLE USE ONLY.

This package contains an HCT/P (Human Cells, Tissues, and Cellular and Tissue Based Products) / CTO (Human Cell, Tissue, Organ) as defined by US Food and Drug Administration (FDA) 21 CFR Part 1271 and the Health Canada Cells, Tissues and Organs for Transplantation Regulations. All tissue has been collected, processed, stored and distributed according to the Standards for Tissue Banking of the American Association of Tissue Banks (AATB), the US FDA regulations, and the Health Canada CTO Regulations and associated Standards (when applicable).

DESCRIPTION / USE

AlloWrap[®] DS and AlloWrap[®] Dry consist of human amniotic membrane that has been processed using a proprietary technology, and is designed with two layers of amniotic tissue with the epithelial layers facing outwards. The products may be used in a variety of procedures as a wound cover or barrier. All tissue is processed and packaged using aseptic technique.

AlloWrap DS is provided in a single patient use package, is moistened with sterile saline solution, and is terminally sterilized using electron beam irradiation.

AlloWrap Dry is provided in a single patient use package, is dehydrated, and is terminally sterilized using electron beam irradiation.

AlloWrap DS and AlloWrap Dry are supplied in multiple sizes for surgical use by licensed clinicians (i.e., physicians, physician's assistants, nurse practitioners).

It is the responsibility of the tissue dispensing service and user (facility/clinician) to maintain the product under appropriate conditions prior to use.

CONTRAINDICATIONS

The presence of infection at the transplantation site is a contraindication for use of tissue allografts.

PRECAUTIONS

- **Restricted to use by a licensed clinician.**
- **Trace amounts of glutaraldehyde may be present and caution should be exercised if the recipient is allergic to this reagent.**
- AlloWrap DS and AlloWrap Dry must be used prior to the expiration date.
- Do not use if the packaging has been damaged and/ or the product has been contaminated. In the event of contamination, discard the product. Damaged packaging should be returned to AlloSource[®].
- Appropriate placement and/or fixation are critical factors in the avoidance of potentially adverse effects.
- As with all biological products, the tissue in AlloWrap DS and AlloWrap Dry has the potential to transmit infectious agents despite processing treatments, extensive donor screening, tissue selection and laboratory tests.
- As with any surgical procedure, the possibility of infection exists.
- Although the production technique is designed to eliminate antigenic properties of the product, the possibility of such a reaction is present.

WARNINGS

- **Human tissue has the potential to transmit infectious agents.** Donor screening, processing treatments and laboratory testing follow stringent specifications to reduce the risk of infectious agent transmission.
- Do not use if the expiration date has been exceeded or if there is evidence of defects in package or label integrity.
- Do not re-sterilize.
- It is the responsibility of the hospital or clinician to maintain tissue for transplantation according to recommended storage conditions. Do not use if tissue has not been stored according to the recommended STORAGE instructions.

ATTENTION:

- Patients receiving any allografts in a surgical procedure should be appropriately informed of the risk associated with these grafts.

POTENTIAL COMPLICATIONS / ADVERSE REACTIONS

Inherent uncertainty exists in medical and social histories and laboratory testing which may not detect known or unknown pathogens. Therefore, the following complications may occur with tissue transplantation:

- Loss of integrity of transplanted tissue with resorption, disintegration, and associated loss of continuity.
- Immune response to transplanted tissue.
- Transmission of known pathogens including Hepatitis B or C, Human T-cell Leukemia/Lymphotropic Virus, Human Immunodeficiency Virus 1 & 2, syphilis or bacteria.
- Transmission or causation of diseases of unknown etiology and characteristics.
- Adverse outcomes potentially attributable to the product must be reported promptly to the manufacturer. Notify AlloSource immediately for return authorization if any dissatisfaction with the product performance or packaging occurs.

DONOR ELIGIBILITY

Donor eligibility (screening and testing) is performed in accordance with US FDA regulations, AATB Standards, and Health Canada CTO regulations. Donor screening includes assessment of the medical and social history as well as physician assessment of the donor to assure that no conditions exist that may make the tissue unacceptable for transplantation. Donor eligibility determination has been made by an AlloSource Medical Director who reviews and approves each donor for processing. Pertinent records may be made available upon request.

SEROLOGICAL TESTING

Communicable disease testing was performed by a laboratory registered with FDA to perform donor testing and certified to perform such testing on human specimens in accordance with the Clinical Laboratory Improvement Amendments of 1988 (CLIA) and 42 CFR Part 493. The testing was conducted using FDA licensed, approved, or cleared donor screening tests. The records of this testing are maintained at AlloSource at the address listed at the bottom of this document. The following required testing was performed and found to be negative or non-reactive.

- Antibody to Human Immunodeficiency Virus 1 & 2 (HIV 1 & 2)
- Human Immunodeficiency Virus Type 1 (HIV-1 NAT)
- Antibody to Hepatitis C (HCV)
- Hepatitis C Virus (HCV NAT)
- Hepatitis B Core IgG/IgM Antibody (HBcAb)
- Hepatitis B Surface Antigen (HBsAg)
- Rapid Plasma Reagin or Serologic Test For Syphilis (RPR or STS)

Additional tests including but not limited to Human T-Cell Lymphotropic Virus Type I & II (HTLV I & II) may have been performed at the time of donor screening and were found to be acceptable for transplantation. A list of additional communicable disease test(s) performed may be provided upon request.

MICROBIAL TESTING

Tissue is subjected to microbiological testing in the course of processing, and must be free of specific aerobic / anaerobic microorganisms and fungal contaminants whose presence would preclude tissue from processing or transplantation.

MEDICAL DIRECTOR ASSESSMENT

Donor eligibility determination is made by the AlloSource Medical Director who reviews and approves each donor for processing. Pertinent records may be made available upon written request.

STORAGE

AlloWrap DS and AlloWrap Dry are packaged either moistened in sterile saline or dehydrated, and should be stored at room temperature in a clean, dry location. **DO NOT REFRIGERATE, FREEZE OR EXPOSE TO EXTREME HEAT.** Tissue shelf life is listed on the product label.

As with all biological materials, some variations in the product should be expected, such as in appearance and in handling.

Do not use any tissue if the expiration date has been exceeded or if there is evidence of defect in the package or label integrity.

DIRECTIONS FOR HANDLING AND PREPARATION OF ALLOWRAP DS

These instructions are intended as guidelines for the use of AlloWrap DS as part of established surgical techniques. They are not intended to replace or change standard procedures or institutional protocols or practices.

All preparation should be performed using aseptic technique.

NOTE: Prepared grafts must be used for the surgical event for which they were intended or otherwise DISCARDED.

1. Graft is contained within 2 pouches. The inner pouch has been sterilized. Using aseptic technique, peel outer pouch and introduce inner-most pouch onto sterile field.
2. Open inner pouch. AlloWrap DS is packaged between two pieces of sterile gauze backing material, moistened with 0.9% sterile saline solution.
3. Wet the graft and backing material with sterile saline solution either in the inner pouch itself or in a sterile preparation container once removed from pouch. If needed, while the tissue is between the two pieces of backing material, cut the graft to desired size using sterile scissors.
4. Using sterile forceps or other surgical tool lift and remove one or both layer(s) of the backing material, then position AlloWrap DS in place. Sterile saline solution may be used to keep the graft moist and pliable. If the backing material is used as a carrier to position the graft into place, remove the remaining layer of the backing material once the graft is secured.

DIRECTIONS FOR HANDLING AND PREPARATION OF ALLOWRAP DRY

These instructions are intended as guidelines for the use of AlloWrap Dry as part of established surgical techniques. They are not intended to replace or change standard procedures or institutional protocols or practices.

All preparation should be performed using aseptic technique.

NOTE: Prepared grafts must be used for the surgical event for which they were intended or otherwise DISCARDED.

1. Graft is contained within 2 pouches. The inner pouch has been sterilized. Using aseptic technique, peel outer pouch and introduce inner-most pouch onto sterile field.
2. Open inner pouch and remove graft from package using sterile forceps or other surgical tool. AlloWrap Dry is packaged dehydrated.
3. If needed, cut the graft to desired size using sterile scissors.
4. Using sterile forceps or other surgical tool, place graft in surgical site or wrap desired tissue with graft. Graft will rehydrate once placed in the patient.

RECORD KEEPING

Regulatory bodies such as the US FDA and Health Canada require that allograft tissue be traceable from the donor to the recipient. The tissue bank (source establishment) is responsible for traceability from the donor to the consignee (transplantation facility, clinician or hospital), and the transplantation facility is responsible for traceability to the recipient. A *Transplantation Record & Feedback Form* and pre-printed peel-off labels are included with each package of tissue. The lot number provided on the labels contains the donor identification number and graft serial number. Record patient identifier, the transplantation facility name and address, the allograft tissue identification information (using the peel-off stickers) and comments regarding the use of the tissue on the *Transplantation Record & Feedback Form*. Return the completed form to AlloSource and retain a copy in the patient medical record. If the tissue has been discarded, please return the *Transplantation Record & Feedback Form* to AlloSource with the graft identification information and reason for discard.

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