

Two layers for double-sided protection

AlloWrap® DS* is a human amniotic membrane designed to provide a biologic barrier following surgical repair. The double-sided membrane features two layers of amniotic tissue orientated with the epithelial layers facing outwards.


ALLOWRAP DS FEATURES

- Packaged wet – no rehydration required
- Non-expansive tissue
- Stable, thin, strong membrane barrier
- Pliable tissue that conforms and settles within the surrounding tissues
- Easy to handle and position into place with or without suturing

CLINICAL BENEFITS

- Immune privileged¹
- Double-sided epithelial layer allows for maximum tissue barrier opportunities as an onlay and/or in wrapping tissue applications
- Tissue remains in surgical site greater than eight weeks²

* AlloWrap DS design is patent pending.



1. Akle CA, Adinolfi M, Welsh KI. Immunogenicity of human amniotic epithelial cells after transplantation into volunteers. Lancet 1981;2:1003-5

2. Samaniego, A.C., Ronholdt, C.J., Ryan, S., Ehrhart, E.J., Turner, S., and Bogdanský, S., Human Amniotic Tissue as an Anti-adhesion, Anti-inflammatory Barrier in an Ovine Spinal Laminectomy Model, AlloSource 2009



*Human amniotic
membrane designed
to provide barrier
protection following
surgical dissection*

AlloWrap® DS

PART CODE	DESCRIPTION
	Sizes (layered double-sided epithelial)
70120004	2 cm x 2 cm
70120008	2 cm x 4 cm
70120016	4 cm x 4 cm
70120032	4 cm x 8 cm

TISSUE PROCESSING

Proprietary aseptic cleaning technology, sterilized and packaged in saline between 2 non-woven protective gauze pads one of which can be used as a carrier

INDICATIONS

Surgical barrier or cover to procedural area used in primary and revision surgical applications for nerves, tendons, vessels and other tissues

SHELF LIFE

Two years, room temperature storage

HANDLING CHARACTERISTICS

Thin, pliable tissue that can be sutured, wrapped or laid within the surgical defect

Ordering Information

TOLL FREE

800. 557. 3587

FAX

720. 873. 0207

AlloSource obtains amniotic tissue for AlloWrap DS only from consenting donors who elect to have scheduled cesarean births. All donors are required to pass stringent medical/social examinations and serologies. All tissue has been collected, processed, stored and distributed according to the Standards for Tissue Banking of the American Association of Tissue Banks (AATB) and Food and Drug Administration (FDA) Regulations.

6278 S Troy Cir
Centennial, CO 80111

MAIN 720. 873. 0213
FAX 720. 873. 0212

allosource.org

