


# AceConnex<sup>®</sup>

## PRE-SUTURED FASCIA

### For Hip Labral Reconstruction and Augmentation

Introducing the next generation labral reconstruction and augmentation allograft, AceConnex, a ready-to-use sterile device, designed specifically for hip arthroscopists to help simplify hip labral procedures.



### SAVES TIME


Pre-suturing saves procedure time by eliminating pre-operative manual allograft preparation in the operating room.



Surgeon time study using fascia and tendons


### FUNCTIONAL

Biomechanical testing shows that AceConnex has similar force distribution as the native labrum.



### CUSTOMIZABLE

Array of pre-sutured fascia sizes, with trimmable regions, that allow for allograft adjustments to match patients' anatomy.



### CONSISTENT

Manufactured pre-sutured fascia ensures consistency, minimizes variability, and remains tightly bound compared to allografts sutured manually in the operating room.\*



SURGEONS HAVE TRUSTED AND IMPLANTED ALLOSOURCE'S FASCIA IN LABRAL RECONSTRUCTION PROCEDURES SINCE 2012

# AceConnex<sup>®</sup>

## PRE-SUTURED FASCIA

### QUALITY\*

Validation testing shows that the device retains functionality at a two year shelf-life.

Each lot is verified with in-process diameter and length measurements.

Safety is ensured through the AlloTrue™ tissue cleansing process and terminal sterilization.



### AceConnex PRE-SUTURED FASCIA

LENGTH	WIDTH	PRESERVATION	REF/PRODUCT #
60 mm	5.0 – 5.4 mm	Frozen	PSF141060-5
60 mm	5.5 – 6.0 mm	Frozen	PSF141060-6
100 mm	5.0 – 5.4 mm	Frozen	PSF141100-5
100 mm	5.5 – 6.0 mm	Frozen	PSF141100-6
140 mm	5.0 – 5.4 mm	Frozen	PSF141140-5
140 mm	5.5 – 6.0 mm	Frozen	PSF141140-6



AceConnex Pre-Sutured Fascia is available in three trimmable lengths and two different diameters, which allows the surgeon to trim the allograft to match the patients' anatomy.

### INDICATIONS FOR USE

The AceConnex Pre-Sutured Fascia is intended for use as a component in soft tissue surgical procedures where constructs, including those with allograft tissue, are used for reconstruction, replacement, or augmentation of the labrum.

\*Data on file

**ORDER NOW**  
800. 557. 3587

AlloSource<sup>®</sup>, one of the largest human tissue providers, honors tissue donors by creating innovative dermis, cartilage, tendon, fascia, bone, and amnion allografts to help heal patients. Since 1994, we have continued to advance our allografts to improve patient outcomes, serving as a trusted tissue partner to the medical community. Learn more at [allosource.org](http://allosource.org).

AceConnex Pre-Sutured Fascia is a 510k cleared medical device. AlloSource is registered with the FDA as a tissue establishment and accredited by the American Association of Tissue Banks.

The AceConnex products are covered under US Patent: 11,607,315

00179-LIT [001] - 06/2024

6278 S Troy Cir  
Centennial, CO 80111  
USA

MAIN 720. 873. 0213

[allosource.org](http://allosource.org)

